a member of the semperit-group

PRODUCT SPECIFICATION

Product name: Spec code: Date of issue:

semperguard® Style NOF-040BK-N-3CZ 28-Aug-2017

Туре	single use examination and disposable protective glove, non sterile	
Labelling	all information on dispenser box	Available sizes
Shape	ambidextrous - straight fingers	
Material	Nitrile Butadiene Rubber (NBR) [not made with natural rubber latex]	S (6-7)
Colour	black	M (7-8)
Inner	online chlorinated / powder free	L (8-9)
Outer	no treatment	XL (9-10)
Cuff/surface	rolled cuff/fingertip textured	
Storage life	3 years	A

Min. thickness [in mm]	Finger	Palm	Cuff
single measured	0.100	0.070	0.050
double measured	0.20	0.14	0.10
	Before aging	After aging	Testing method
Elongation	min. 500 %	min. 400 %	ASTM D6319
Tensile strength	min. 14 MPa	min. 14 MPa	ASTM D6319
AQL (for pinholes)	≤ 1.5		EN 374-2 EN 455-1
Residual powder / Powder content	≤ 2 mg/glove		EN 455-3
Force at break	median ≥ 6.0 N		EN 455-2
Length	median 240 mm		EN 455-2

Storage instructions

Store according to ISO 2230:2002. Ensure that the area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper-ions discolour the glove. Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.

> semperguard® Style NOF-040BK-N-3CZ_Version D_August 2017_EN Replaces all previous versions and valid for 2 years after date of issue

Page 1/2

Standards, guidelines & quality certificates

Quality certification	ISO 9001 ISO 13485 ISO 14001 OSHAS 18001
Conformity to directives	Personal Protective Equipment 89/686/EEC: Cat III for complex risks
	Medical Device Directive 93/42/EEC: Class I
Confomity to standards	EN 420 EN 455 1-4 EN 374-1-3
	ASTM D6319 ASTM F1671
Conformity to food regulation	1935/2004 (EC)
labelling according to	ISO 15223 EN 420 EN 374

Warning and reporting system

Potentially allergenic ingredients

This product contains accelerators (Dithiocarbamate type, Zinc-mercaptobenzothiazol) not to be used in a hypersensitivity of these substances.

For further information, the actual version of the ingredients list is given on request.

Medical device vigilance and reporting system

According to the official reporting criteria of MDD, incidents caused by examination gloves or surgical gloves must be reported immediately to our medical device reporting officer. E-Mail: sempermed.complaints@semperitgroup.com or Tel.: +43 2630 310 ext. -510; Fax: -549

Authorized signator

C. Eichler

Page 2/2

Product management

A. Scharl

Regulatory affairs

J. Glantschnig

semperguard® Style NOF-040BK-N-3CZ_Version D_August 2017_EN Replaces all previous versions and valid for 2 years after date of issue