

RefectoCil tint remover 150 ml

Description

» Color stain remover

Removes stains that develop during tinting easily and quickly.

Usage:

Immediately after tinting, dampen a cotton ball or a Q-tip with RefectoCil Tint remover and gently rub off the stain.

The quicker the stain is treated the easier its removal.

Tip!

The quicker the stain is treated the easier its removal. Only suited for RefectoCil oxidative tints! (article 6160-6184)

Ingredients:

Aqua, Isopropyl Alcohol, Polysorbate 80, Citrus Limonum, Limonene, Litsea Cubeba Fruit Oil, Citral, Linalool, Citronellol, Geraniol.

For more information about the assortment and tutorials, please visit Refecto Cil's Youtube Channel [here](#)